

Backing-up and Restoring 6.x Office Medic Database

This hot tip assumes you have installed Office Medic into its default directory.

Database Back-up Instructions

Backing-up your database protects you from losing your patient data should a catastrophic event occur. Regular back-ups of the database should be maintained. Follow the steps below to back-up the database:

1. Close Office Medic.
2. Open folder: C:\Vectraplex\Database.
3. Copy the two files VectraplexECG.MDF and VectraplexECG_Log.LDF to a secure location. This is the back-up copy of your Office Medic database. Copy these files as often as needed to maintain a current back-up file.

Database Restore Instructions

Warning! Restoring the database erases all of the data located in Office Medic and replaces it with the data contained in the back-up file. Data that was acquired after the date of the last back-up will be lost and cannot be recovered.

Follow the steps below to restore the database:

1. Close Office Medic.
2. Copy and paste the two back-up files into the following location: C:\Vectraplex\Database.
3. Open Office Medic.

The database should look exactly as it did on the date of the last back-up.