

EasyABPM software can be used for programming, analyzing, interpreting and reporting of data recorded by ABPM-05. The said software is the right choice for active **general practitioners**, providing:

1. **Automated Statistics for all Periods**

EasyABPM software solution displays statistics for all (active, passive, morning and special) periods. Statistics include systolic blood pressure, diastolic blood pressure, mean arterial pressure (MAP), pulse pressure (PP), pulse rate, minimum and maximum values, standard deviation (SD), percent time elevation (PTE), hypertension load (hyperbaric index), diurnal index (day/night index) and morning surge.

2. **Default Measurement Plans**

Loadable default ABPM measurement plans make it easy to configure any Meditech 24 hour ABPM quickly.

3. **Customizable Sampling Periods**

4 individual sampling periods (active or awake, passive or asleep, morning and special) are available.

4. **Graphic and Tabular Presentation of 24 hour ABPM Data**

A visual blood pressure chart with blood pressure plotted on the vertical axis, time on the horizontal axis with systolic, diastolic values and pulse are available in the software. All data are displayed in tabular or chart format as well.

5. **Hourly Averages**

Hourly blood pressure averages are displayed to support patient treatment.

6. **Editable Patient and Study Info**

Freely editable patient and examination information with the possibility of writing notes and opinions also, to show physical conditions or patient history in the report.

7. **Integrated Database**

Easy-to-use, integrated patient database.

8. **Data List**

Interpretive 24-hour ABPM data are recorded automatically.

9. **Data Limits**

The user can freely define blood pressure threshold limits for statistical calculations.

10. **Printable Reports**

Multiple report configurations, editable reports or short, standard and full versions are available.