

# Q200/HE System Requirements

Currently Holter LX Analysis version 5.4CQ ships with the Q200/HE. See below for system requirements for the Q200/HE.

## Minimum System Requirements for Desktop, Laptop and Tablet PCs:

<b>Operating System:</b>	Microsoft® Windows® XP or Windows 7
<b>Free Disk Space:</b>	10 GB
<b>RAM:</b>	1 GB
<b>Processor:</b>	1 GHz or faster
<b>Screen Resolution:</b>	1024 x 768
<b>Connectivity:</b>	USB flashcard reader (included) or PC Card Slot